

FILE COPY

4

AD

TECHNICAL REPORT ARCCB-TR-88033

**DETERMINATION OF ETHYLENE GLYCOL DEGRADATION
PRODUCTS IN CHROMIUM PLATING AND ASSOCIATED
POLISHING SOLUTIONS BY ION CHROMATOGRAPHY**

AD-A200 185

SAMUEL SOPOK

AUGUST 1988

DTIC
ELECTE
OCT 1 2 1988
S H



**US ARMY ARMAMENT RESEARCH,
DEVELOPMENT AND ENGINEERING CENTER
CLOSE COMBAT ARMAMENTS CENTER
BENÉT LABORATORIES
WATERVLIET, N.Y. 12189-4050**



APPROVED FOR PUBLIC RELEASE; DISTRIBUTION UNLIMITED

88 1011 212

REPORT DOCUMENTATION PAGE		READ INSTRUCTIONS BEFORE COMPLETING FORM
1. REPORT NUMBER ARCCB-TR-88033	2. GOVT ACCESSION NO.	3. RECIPIENT'S CATALOG NUMBER
4. TITLE (and Subtitle) DETERMINATION OF ETHYLENE GLYCOL DEGRADATION PRODUCTS IN CHROMIUM PLATING AND ASSOCIATED POLISHING SOLUTIONS BY ION CHROMATOGRAPHY		5. TYPE OF REPORT & PERIOD COVERED Final
		6. PERFORMING ORG. REPORT NUMBER
7. AUTHOR(s) Samuel Sopok		8. CONTRACT OR GRANT NUMBER(s)
9. PERFORMING ORGANIZATION NAME AND ADDRESS US Army ARDEC Benet Laboratories, SMCAR-CCB-TL Watervliet, NY 12189-4050		10. PROGRAM ELEMENT, PROJECT, TASK AREA & WORK UNIT NUMBERS AMCMS No. 6126.23.1BL0.0 PRON No. 1A72ZJ7BNM\$C
11. CONTROLLING OFFICE NAME AND ADDRESS US Army ARDEC Close Combat Armaments Center Picatinny Arsenal, NJ 07806-5000		12. REPORT DATE August 1988
		13. NUMBER OF PAGES 5
14. MONITORING AGENCY NAME & ADDRESS (if different from Controlling Office)		15. SECURITY CLASS. (of this report) UNCLASSIFIED
		15a. DECLASSIFICATION/DOWNGRADING SCHEDULE
16. DISTRIBUTION STATEMENT (of this Report) Approved for public release; distribution unlimited.		
17. DISTRIBUTION STATEMENT (of the abstract entered in Block 20, if different from Report)		
18. SUPPLEMENTARY NOTES Submitted to <u>LC-GC Magazine of Liquid and Gas Chromatography</u> .		
19. KEY WORDS (Continue on reverse side if necessary and identify by block number) Ethylene Glycol Acid Degradation Products Chromium Plating Solutions Polishing Solutions Ion Chromatography		
20. ABSTRACT (Continue on reverse side if necessary and identify by block number) Ethylene glycol resulting from cooling system leaks can adversely affect plating properties when added to chromium plating and associated polishing solutions. Ion chromatography can be used to monitor these leaks by quantitatively determining the glycolic, oxalic, and formic acid degradation products of ethylene glycol.		

TABLE OF CONTENTS

	<u>Page</u>
INTRODUCTION	1
EXPERIMENTAL PROCEDURE	1
RESULTS AND DISCUSSION	2
REFERENCES	4

TABLES

I. DETERMINATION OF DEGRADATION PRODUCTS OF ETHYLENE GLYCOL	3
---	---



Accession For		<input checked="" type="checkbox"/>
NTIS GRA&I		
DTIC TAB		
Unannounced Justification		
By _____		
Distribution/		
Availability		
Special		
Dist		
A-1		

INTRODUCTION

Many general methods (refs 1,2) exist for determining organic acids by high performance ion chromatography exclusion (HPICE). Past investigations (refs 3,4) have shown that ion chromatography exclusion is useful for determining the organic acid degradation products of ethylene glycol used in solar energy collectors. The oxidative properties of both chromium plating and associated polishing solutions also produce organic acid degradation products of ethylene glycol which compromise plating quality (refs 5,6).

A method is given here using HPICE for monitoring and determining organic acid degradation products due to ethylene glycol coolant leaks in chromium plating and associated polishing solutions.

EXPERIMENTAL PROCEDURE

The analysis system is the Dionex model 2020i ion chromatograph and the Dionex model 4270 integrator (attenuation is 1024). The system flowstream consists of the eluent reservoir (0.001 M hydrochloric acid), eluent delivery system (1.0 ml/min), injector valve (100 μ l), separator column (Dionex, HPICE-AG1), suppressor column (Dionex, ISC), and conductivity detector (Dionex, 10 μ S).

Calibration standards are composed of only glycolic, oxalic, and formic acids since the presence of sulfuric, phosphoric, or chromic acids would decompose the organic acids.

Samples are composed of 25 ml of reagent grade ethylene glycol and 975 ml of either chromium plating or associated polishing solutions. The chromium plating solution contains 2.50 g/l sulfuric acid and 250 g/l chromic acid, while

References are listed at the end of this report.

the polishing solution contains 680 g/l phosphoric acid and 840 g/l sulfuric acid. The sample solutions are each allowed to react with ethylene glycol for two hours under stirred conditions allowing gases to evolve. Although vigorous evolution of carbon dioxide subsides after ten minutes, these samples are stirred for two hours to assure that degradation products equilibrate in these metal finishing solutions. After stirring, each sample requires a 250 dilution before analysis.

RESULTS AND DISCUSSION

In Table I, the first three solutions provide calibration data for glycolic, oxalic, and formic acids which are organic acid degradation products of ethylene glycol in chromium plating and associated polishing solutions. Solutions four and five in Table I contain the respective results of chromium plating and associated polishing solutions for determinations of these organic acid products of ethylene glycol. For all solutions in Table I, the retention times of glycolic, oxalic, and formic acids are 11.7, 7.0, and 12.8 minutes, respectively.

TABLE I. DETERMINATION OF DEGRADATION PRODUCTS OF ETHYLENE GLYCOL

Solution	Glycolic Acid		Oxalic Acid		Formic Acid	
	Conc.(ppm)	H(μ S)	Conc.(ppm)	H(μ S)	Conc.(ppm)	H(μ S)
1	10	2.61	10	4.09	10	3.53
2	20	5.19	20	8.19	20	7.00
3	30	7.71	30	11.93	30	10.50
4	2	0.40	12	4.85	15	5.41
5	4	1.11	19	7.62	26	9.10

As suggested elsewhere (refs 3-5), the oxidative degradation of ethylene glycol has ordered steps in these metal finishing solutions that include glycolic, oxalic, formic, and carbonic acids. In these chromium plating solutions, trivalent chromium is a by-product of the oxidation of ethylene glycol. The composition of ethylene glycol degradation products depends on acid composition, acid concentration, reaction time, and reaction temperature.

A contributing factor in the success of this method is that HPICE does not retain cations, sulfate, phosphate, or chromate, and no interference exists between these non-retained species and the retained organic acid analytes.

In addition to the analytes mentioned, there is a very small peak at 15.3 minutes with a height of 0.1 μ S that appears to be carbonic acid, but this is difficult to investigate in acidic solutions and eluent. There is a great deal of carbonic acid liberated by the oxidation process and it is reasonable to assume that some is still present under weak acidic conditions.

This method provides a baseline resolution of analyte peaks and quantitatively determines these analytes to parts per million levels in these metal finishing solutions.

REFERENCES

1. F. Smith and R. Chang, The Practice of Ion Chromatography, John Wiley and Sons, New York, 1983.
2. J. Fritz, D. Gjerde, and C. Pohlandt, Ion Chromatography, Huethig, Heidelberg, Germany, 1982.
3. W. Rossiter, M. Godette, P. Brown, and K. Galuk, "An Investigation of the Degradation of Aqueous Ethylene Glycol and Propylene Glycol Solutions Using Ion Chromatography," Solar Energy Materials, Vol. 11, 1985, pp. 455-467.
4. W. Rossiter, P. Brown, and M. Godette, "The Determination of Acidic Degradation Products in Aqueous Ethylene Glycol and Propylene Glycol Solutions Using Ion Chromatography," Solar Energy Materials, Vol. 9, 1983, pp. 267-279.
5. R. Morrison and R. Boyd, Organic Chemistry, Third Edition, Allyn and Bacon, Inc., Boston, 1974, pp. 528, 634.
6. "Metal Finishing Handbook," Metals and Plastics Publication, Inc., Hackensack, NJ, 1982.

TECHNICAL REPORT INTERNAL DISTRIBUTION LIST

	NO. OF COPIES
CHIEF, DEVELOPMENT ENGINEERING BRANCH	
ATTN: SMCAR-CCB-D	1
-DA	1
-DC	1
-DM	1
-DP	1
-DR	1
-DS (SYSTEMS)	1
CHIEF, ENGINEERING SUPPORT BRANCH	
ATTN: SMCAR-CCB-S	1
-SE	1
CHIEF, RESEARCH BRANCH	
ATTN: SMCAR-CCB-R	2
-R (ELLEN FOGARTY)	1
-RA	1
-RM	1
-RP	1
-RT	1
TECHNICAL LIBRARY	5
ATTN: SMCAR-CCB-TL	
TECHNICAL PUBLICATIONS & EDITING UNIT	2
ATTN: SMCAR-CCB-TL	
DIRECTOR, OPERATIONS DIRECTORATE	1
ATTN: SMCWV-OD	
DIRECTOR, PROCUREMENT DIRECTORATE	1
ATTN: SMCWV-PP	
DIRECTOR, PRODUCT ASSURANCE DIRECTORATE	1
ATTN: SMCWV-QA	

NOTE: PLEASE NOTIFY DIRECTOR, BENET LABORATORIES, ATTN: SMCAR-CCB-TL, OF ANY ADDRESS CHANGES.

TECHNICAL REPORT EXTERNAL DISTRIBUTION LIST

	<u>NO. OF COPIES</u>		<u>NO. OF COPIES</u>
ASST SEC OF THE ARMY RESEARCH AND DEVELOPMENT ATTN: DEPT FOR SCI AND TECH THE PENTAGON WASHINGTON, D.C. 20310-0103	1	COMMANDER ROCK ISLAND ARSENAL ATTN: SMCRI-ENM ROCK ISLAND, IL 61299-5000	1
ADMINISTRATOR DEFENSE TECHNICAL INFO CENTER ATTN: DTIC-FDAC CAMERON STATION ALEXANDRIA, VA 22304-6145	12	DIRECTOR US ARMY INDUSTRIAL BASE ENGR ACTV ATTN: AMXIB-P ROCK ISLAND, IL 61299-7260	1
COMMANDER US ARMY ARDEC ATTN: SMCAR-AEE	1	COMMANDER US ARMY TANK-AUTMV R&D COMMAND ATTN: AMSTA-DDL (TECH LIB) WARREN, MI 48397-5000	1
SMCAR-AES, BLDG. 321	1	COMMANDER US MILITARY ACADEMY ATTN: DEPARTMENT OF MECHANICS WEST POINT, NY 10996-1792	1
SMCAR-AET-O, BLDG. 351N	1		
SMCAR-CC	1		
SMCAR-CCP-A	1		
SMCAR-FSA	1		
SMCAR-FSM-E	1	US ARMY MISSILE COMMAND REDSTONE SCIENTIFIC INFO CTR ATTN: DOCUMENTS SECT, BLDG. 4484 REDSTONE ARSENAL, AL 35898-5241	2
SMCAR-FSS-D, BLDG. 94	1		
SMCAR-IMI-I (STINFO) BLDG. 59	2		
PICATINNY ARSENAL, NJ 07806-5000			
DIRECTOR US ARMY BALLISTIC RESEARCH LABORATORY ATTN: SLCBR-DD-T, BLDG. 305 ABERDEEN PROVING GROUND, MD 21005-5066	1	COMMANDER US ARMY FGN SCIENCE AND TECH CTR ATTN: DRXST-SD 220 7TH STREET, N.E. CHARLOTTESVILLE, VA 22901	1
DIRECTOR US ARMY MATERIEL SYSTEMS ANALYSIS ACTV ATTN: AMXSY-MP ABERDEEN PROVING GROUND, MD 21005-5071	1	COMMANDER US ARMY LABCOM MATERIALS TECHNOLOGY LAB ATTN: SLCMT-IML (TECH LIB) WATERTOWN, MA 02172-0001	2
COMMANDER HQ, AMCCOM ATTN: AMSMC-IMP-L ROCK ISLAND, IL 61299-6000	1		

NOTE: PLEASE NOTIFY COMMANDER, ARMAMENT RESEARCH, DEVELOPMENT, AND ENGINEERING CENTER, US ARMY AMCCOM, ATTN: BENET LABORATORIES, SMCAR-CCB-TL, WATERVLIET, NY 12189-4050, OF ANY ADDRESS CHANGES.

TECHNICAL REPORT EXTERNAL DISTRIBUTION LIST (CONT'D)

	<u>NO. OF COPIES</u>		<u>NO. OF COPIES</u>
COMMANDER US ARMY LABCOM, ISA ATTN: SLCIS-IM-TL 2800 POWDER MILL ROAD ADELPHI, MD 20783-1145	1	COMMANDER AIR FORCE ARMAMENT LABORATORY ATTN: AFATL/MN EGLIN AFB, FL 32542-5434	1
COMMANDER US ARMY RESEARCH OFFICE ATTN: CHIEF, IPO P.O. BOX 12211 RESEARCH TRIANGLE PARK, NC 27709-2211	1	COMMANDER AIR FORCE ARMAMENT LABORATORY ATTN: AFATL/MNF EGLIN AFB, FL 32542-5434	1
DIRECTOR US NAVAL RESEARCH LAB ATTN: MATERIALS SCI & TECH DIVISION CODE 26-27 (DOC LIB) WASHINGTON, D.C. 20375	1 1	METALS AND CERAMICS INFO CTR BATTELLE COLUMBUS DIVISION 505 KING AVENUE COLUMBUS, OH 43201-2693	1

NOTE: PLEASE NOTIFY COMMANDER, ARMAMENT RESEARCH, DEVELOPMENT, AND ENGINEERING CENTER, US ARMY AMCCOM, ATTN: BENET LABORATORIES, SMCAR-CCB-TL, WATERVLIET, NY 12189-4050, OF ANY ADDRESS CHANGES.